

Final Project Presentation

EMF Wearable

A Solution for Children with Neuroblastoma

*By Khanh Nguyen, Manish Sharma, & Zachary Kunz*

# Outline

## I. Project Background

## II. EMF Wearable: A Solution for Neuroblastoma

## III. EMF Wearable

A. Circuit Design with EAGLE


B. Device Design with Fusion 360

C. Device Simulation with COMSOL Multiphysics & Sim4Life

D. Safety issues

## IV. IP & FDA Issues


## V. Future Development



[15]

# I. Project Background


- ***NEED STATEMENT: A diagnosis and treatment tool for neuroblastoma to reduce average hospitalization stay time and cost.***
- Neuroblastomas are one of the most prevalent types of pediatric cancer
- Typically latches around the adrenal medulla or paraspinal sites
- Current treatment options include
  - Chemotherapy - in more advanced cases
    - High-dose myeloablative chemotherapy
  - Surgery
    - Autologous hematopoietic stem cell transplantation
  - Medical devices



# I. Project Background

- Average total cost incurred by families suffering from Neuroblastoma over the course of treatment was **\$205,747**
- Since the year 2019 the incidence of neuroblastoma occurring has increased dramatically
  - Drug market size will grow by 190.82 million US
- Emotional and Financial burden from pediatric cancer
  - Pushes around 15% of families experiencing this under the poverty line
  - Loss of jobs due to extensive traveling and caregiving times

## II. EMF Wearable for Neuroblastoma




**Friendly electrode design**

**One device, two modes of operations:** ROS at low frequencies & TTF at high frequencies

**Programmable & customizable** (based on disease conditions)

**Safe & Effective:** High efficiency & safe for use on children and at home

## II. EMF Wearable for Neuroblastoma



**Mode 1: 200-300Hz EMF**


Tumor deep inside the body

**Mode 2: 200kHz-300kHz EMF**

Tumor near the skin or in sensitive tissue

Combination of mode 1 & 2 is also possible

# III. EMF Wearable - Circuit Design with EAGLE



The circuit consists of:

- Coil
- Resistor
- Capacitor
- AC Power Supply

### III. EMF Wearable - Device Design with Fusion 360




#### Coil Dimensions

- ◆ 5 revolution coil
- ◆ Outer radius of 30 mm
- ◆ Inner radius is 4 mm
- ◆ Wire radius 1 mm



#### Box Dimensions

- ◆ 70 mm x 70 mm x 18 mm
- ◆ 25 holes with 0.5 mm radii






### III. EMF Wearable - Device Design with Fusion 360

#### Materials


- ◆ Enclosure - Soft polymer
- ◆ Coil - Copper
- ◆ Wires - Plastic coated insulated wires

#### Layout

- ◆ 5 mm on all sides away from coil to prevent burns
- ◆ 25 holes to allow the heat to escape
- ◆ 4 prongs on top of box to keep lid in place




### III. EMF Wearable - Simulation with COMSOL Multiphysics



- Interested Results
  - Magnetic Field Strength & Flux Density
  - Electric Field Strength
  - Thermal Distribution
- Model Parameters
  - Coil shape and size
  - Enclosure materials
  - Flowing current & applied power
  - Applied frequency


### III. EMF Wearable - Simulation with COMSOL Multiphysics






- EMF components vary with
  - Coil shape and size
  - Enclosure materials
  - Flowing current & applied power
  - Applied frequency
- COMSOL Multiphysics (Ver 5.2) Limitations
  - Unavailability of RF module

=> Can't simulate at different frequencies

=> *Sim4Life is used for this purpose instead*




# III. EMF Wearable - Simulation with Sim4Life




- Sim4Life is a software that specializes in physics simulations with biological models
- Sim4Life is especially useful with simulation of low frequencies (under 1GHz)
- Since the RF module is NOT available in COMSOL Multiphysics ver 5.2, Sim4Life is an alternative. However, some of the main disadvantages for this option is
  - Can not simulate model with high frequencies
  - The student version (Sim4Life Light) does not support 3D view of the analysis
  - Limiting computation capacity compared to COMSOL Multiphysics

### III. EMF Wearable - Simulation with Sim4Life




- Model Parameters
  - Coil shape and size
    - Shape: Circular, spiral
    - Inner radius: 10mm
    - Outer Radius: 30mm
    - Turns: 5
  - Enclosure materials: Soft polymer
  - Applied Current: 0.5A
  - Applied frequency: 200kHz


# III. EMF Wearable - Simulation with Sim4Life




Magnetic field strength




Magnetic flux density



Electric field strength






Magnetic field strength - Multislices



Magnetic field strength - Vector View


Every component of the EMF can be analyzed  
Data is shown mostly in slice and surface views

### III. EMF Wearable - Simulation with Sim4Life




- Data can be extracted & plotted along a 2D line
- Detail information regarding the field strength and its location can be extracted from the plot (in .csv file)


# III. EMF Wearable - Simulation with Sim4Life






Magnetic field strength



Magnetic flux density



Electric field strength



At 2.2cm above the center of the coil

- Magnetic field strength:  
15.07A/m
- Magnetic flux density:  
 $1.89 \times 10^{-5}$  T
- Electric field strength:  
 $1.89 \times 10^{-8}$  V/m

Data extracted


Those results can be easily simulated in COMSOL Multiphysics using “derived values”



### III. EMF Wearable - Safety Issues

#### *Thermal Dissipation*

- Holes on the coil chamber to provide an inward flow of air and gives natural cooling
- Vapor cooling: Vapor cooling is employed when high power and heat fluxes are given to a small evaporator area. The heat input into the evaporator vaporizes liquid within the evaporator wick during operation. The vapor then circulates around the chamber forming an isothermal heat spreader. The vapor condenses on the condenser surfaces, where heat is removed by the convection of natural cooling.



[16]

#### *Electrical Safety*

- There is a switch with a thermal sensor (which can be calibrated), so the button will auto cut the power supply when the temperature rises above a specific limit.

# III. EMF Wearable - Safety Issues

## EMF Exposure Potential Harmful Effects

- Possible tissue damage due to overheating at high frequencies (>100kHz)
- Possible nerve stimulation & cell membrane permeability change at low frequencies (<100kHz)
- Carcinogenic? - Not enough evidence in humans
- Interference with electronic medical devices
- Safety guidelines have been setup to avoid potential harmful effects of EMF
  - **ICNIRP 2010** by the International Commission on Non-Ionizing Radiation Protection
  - **IEEE C95.1** by the International Committee on Electromagnetic Safety

GENERAL PUBLIC EXPOSURE LIMITS <sup>c</sup>			
Exposure of	E-field strength $E$ (V / m)	Magnetic field strength $H$ (A / m)	Magnetic flux density $B$ (μT)
head and torso		163	205
limbs		900	1130
whole body	614		
CONTROLLED ENVIRONMENT EXPOSURE LIMITS <sup>d</sup>			
Exposure of	E-field strength $E$ (V / m)	Magnetic field strength $H$ (A / m)	Magnetic flux density $B$ (μT)
head and torso		490	615
limbs		900	1130
whole body	1842		

EMF Exposure Limit by IEEE C95.1 [14]



### Basic restrictions

Parameter	Frequency range	ΔT	Spatial averaging	Temporal averaging	Health effect level	Reduction factor	Workers	Reduction factor	General public
Core ΔT	100 kHz-300 GHz	1°C	WBA*	30 min	4 W/kg	10	0.4 W/kg	50	0.08 W/kg
Local ΔT (Head & Torso)	100 kHz-6 GHz	2°C	10 g	6 min	20 W/kg	2	10 W/kg	10	2 W/kg
Local ΔT (Limbs)	100 kHz-6 GHz	5°C	10 g	6 min	40 W/kg	2	20 W/kg	10	4 W/kg
Local ΔT (Head & Torso, Limbs)	>6-300 GHz 30-300 GHz	5°C	4 cm <sup>2</sup> 1 cm <sup>2</sup>	6 min 6 min	200 W/m <sup>2</sup> 400 W/m <sup>2</sup>	2	100 W/m <sup>2</sup> 200 W/m <sup>2</sup>	10	20 W/m <sup>2</sup> 40 W/m <sup>2</sup>

\* WBA: whole body average

SAR Limit by ICNIRP 2020 [13]

GENERAL PUBLIC EXPOSURE LIMITS <sup>a</sup>			
Frequency range	E-field strength $E$ (V / m)	Magnetic field strength $H$ (A / m)	Magnetic flux density $B$ (μT)
3 kHz – 10 MHz	83	21	27
OCCUPATIONAL EXPOSURE LIMITS <sup>b</sup>			
Frequency range	E-field strength $E$ (V / m)	Magnetic field strength $H$ (A / m)	Magnetic flux density $B$ (μT)
3 kHz – 10 MHz	170	80	100

EMF Exposure Limit by ICNIRP 2010 [14]

## IV. IP Issues

- Utility Patent - A patent that covers the creation of a new or improved product and prohibits others from using or selling the product
  - Lasts 20 years
  - Claims - The usage of an electromagnetic field for targeted therapy of cancer cells by reactive oxygen species activation using frequencies 200-300Hz in combination with tumor treating field of 200-300kHz
- A similar device called Optune Wearable (NovoFFT-100A)
  - Has the classifications A61N1/32, A61K31/196, A61N1/36002, A61N1/40, and A61P43/00
- Freedom to operate: If we have the ability to have very different technology to optune than we will be able to patent but if our technology is similar we will not be able to patent.
  - Our idea has novelty as compared to Optunes design due to the fact that we have a different design and technical detail from them. This would allow for our company to have freedom to operate without paying royalties to Optune

# IV. FDA Regulation

EMF Wearable is a class III device and will go through Pre-Market Approval

Clinical studies will be conducted with special cautions regarding use at home and on children

## Clinical Study Design Details

- Disease: Pathological or histological evidence of Neuroblastoma
- Age range: children aged 1-5 years old
- Population: 1000 subject of diverse ethnic, racial, and geographical background
- Cancer stage: Stage I-IV
- Other disease diagnosis/treatment specifications:
  - Newly diagnosed or recurrent
  - Treated in conjunction with other therapies (if any)
- Life expectancy: At least 3 months from the beginning of the study
- Informed consent: Parents' consent is mandatory
- Access to the device: Having a parent or a caregiver being able to operate the device independently
- Testing environment: Uncontrolled environment (for at-home usage)
- Result: Cancer progress after treatment

# V. Future Development

## ***Device Design***

Incorporate array design from single-coil to customize the EMF exposure (shape, size, and exposure level)

Improve thermal dissipation efficiency

Create a software that assist the users and report data to physicians

## ***Device Testing***

Complete simulation in COMSOL Multiphysics & Sim4Life

Evaluate the performance of the device in real-life setting

Evaluate the overall safety of the device considering general guidelines regarding EMF exposure



# References

1. “How Optune® Works,” *Optune.com*. [Online]. Available: <https://www.optune.com/how-optune-works>. [Accessed: 06-Mar-2022].
2. “Novocure’s Tumor Treating Fields: Innovative brain cancer therapy with survival and safety benefits,” *Nature news*. [Online]. Available: <https://www.nature.com/articles/d42473-018-00156-3>. [Accessed: 06-Mar-2022].
3. Sengupta, & Balla, V. K. (2018). A review on the use of magnetic fields and ultrasound for non-invasive cancer treatment. *Journal of Advanced Research*, 14, 97–111. <https://doi.org/10.1016/j.jare.2018.06.003>
4. S. Tsubota and K. Kadomatsu, “Origin and initiation mechanisms of neuroblastoma,” *Cell and Tissue Research*, vol. 372, no. 2, pp. 211–221, 2018.
5. Cheung NK, Dyer MA (2013) Neuroblastoma: developmental biology, cancer genomics and immunotherapy. *Nat Rev Cancer* 13:397–411
6. A. B. Mariotto, K. Robin Yabroff, Y. Shao, E. J. Feuer, and M. L. Brown, “Projections of the cost of cancer care in the United States: 2010-2020,” *JNCI Journal of the National Cancer Institute*, vol. 103, no. 2, pp. 117–128, 2011.
7. K. Bona, V. Dussel, L. Orellana, T. Kang, R. Geyer, C. Feudtner, and J. Wolfe, “Economic impact of advanced pediatric cancer on families,” *Journal of Pain and Symptom Management*, vol. 47, no. 3, pp. 594–603, 2014.
8. M. Bilodeau, C. Ma, H. Al-Sayegh, J. Wolfe, and K. Bona, “Household material hardship in families of children post-chemotherapy,” *Pediatric Blood & Cancer*, vol. 65, no. 1, 2017.
9. K. Roser, F. Erdmann, G. Michel, J. F. Winther, and L. Mader, “The impact of childhood cancer on parents' socio-economic situation—a systematic review,” *Psycho-Oncology*, vol. 28, no. 6, pp. 1207–1226, 2019.
10. *Management*, vol. 47, no. 3, pp. 594–603, 2014.  
FDA. Center for Devices and Radiological Health, “Overview of device regulation,” *U.S. Food and Drug Administration*. [Online]. Available: <https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/overview-device-regulation>. [Accessed: 13-Mar-2022].
11. FDA, “Premarket approval (PMA),” *accessdata.fda.gov*. [Online]. Available: <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpma/pma.cfm?id=P100034S013>. [Accessed: 13-Mar-2022].
12. FDA, “Humanitarian device exemption (HDE),” *accessdata.fda.gov*. [Online]. Available: <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfhde/hde.cfm?id=H180002>. [Accessed: 13-Mar-2022].
13. ICNIRP, “Guidelines for limiting exposure to electromagnetic fields (100 KHz to 300 GHz),” *Health Physics*, vol. 118, no. 5, pp. 483–524, 2020.
14. T. M. Fisher, K. B. Farley, Y. Gao, H. Bai, and Z. T. H. Tse, “Electric vehicle wireless charging technology: a state-of-the-art review of magnetic coupling systems,” *Wireless Power Transfer*, vol. 1, no. 2, pp. 87–96, 2014.
15. “Disposable foam electrodes” *Gima, Articoli Medicali, apparecchi per Medicina, Elettromedicali, endoscopi, Misuratori Medicali, Arredo Ambulatorio, Elettrobisturi, articoli per il Pronto Soccorso, articoli per la sterilizzazione*. [Online]. Available: [https://www.gimaitaly.com/prodotti.asp?sku=33373&dept\\_selected=58&dept\\_id=585](https://www.gimaitaly.com/prodotti.asp?sku=33373&dept_selected=58&dept_id=585). [Accessed: 01-May-2022].
16. *Vapor Chambers and their use in thermal management (part 2 of 2)*. Vapor Chambers and Their Use in Thermal Management (part 2 of 2). (n.d.). Retrieved May 2, 2022, from <http://qats.com/cms/2010/12/10/vapor-chambers-and-their-use-in-thermal-management-part-2-of-2/>